

Aminoallyl-dUTP

Aminoallyl-dUTP

Absorption max.:
260 nm (in water)

Molar absorbance:
9.000 M⁻¹cm⁻¹


Comments:

- Concentration: 1 mg/ml
- Suitable for microarray experiments, FisH microscopy

Available Modifications

Available Modifications	Molecular Weight (g • mol ⁻¹)	Molecular formula	Productnumber
Aminoallyl-dUTP	547.66	C ₁₂ H ₁₆ N ₃ O ₁₄ P ₃ * 4Li	010-00

Formel



Order

Available Modifications	Amount	Productnumber	Price	Enter quantity
Aminoallyl-dUTP	1 mg	010-00	100,00 €	